

AEROCASE® - RPL1C

Emergency backpack 1R in "Large" made of AEROTex® - PLAN (Tarpaulin) or AEROTex® - PES (Polyester).

- Main compartment with stable body frame and removable partition
- 5 transparent, dimensionally stable modular bags with hook and loop fastening in the interior: 3 in the main compartment, 2 in the unfoldable front
- Elastic fixation straps inside, 2 side pockets, 2 front pockets, transparent compartment
- Ergonomically improved carrying handle and padded back carrying straps with quick release straps for fixation to the body
- Reinforced bottom protection made of tarpaulin material
- Bidirectional zippers
- Hook and loop strip for patches



Padded, ergonomic back carrying straps

Mesh in the back area



For all models



High quality black AEROTex® - PES (Polyester) in the interior of the main compartment and the front pocket

Particularly tear-resistant and easily washable! - for all models -



Clear interior concept

Individual organization

- 5 modular bags in different sizes
- Tarpaulin - single colour (Red / Blue / Luminous Yellow)
- Polyester - colour-coded
- Transparent window and handles



Model: Tarpaulin




Model: Polyester

AERotex® - PLAN



Material	Colour	Size	Dimensions (W x H x D)	Weight	Equipment	Item No.
AERotex® - PLAN (Tarpaulin)	■ Signal Red (RAL 3001)	L	approx. 36 x 47 x 26 cm	approx. 3 kg	Modular Bag Set, Tarpaulin, Red (Tarpaulin)	① HT03-RPL1C-R
	■ Luminous Yellow (RAL 1026)				Modular Bag Set, Tarpaulin, Blue	② HT03-RPL1C-B
	■ Ultramarine Blue (RAL 5002)				Modular Bag Set, Tarpaulin, Luminous Yellow	③ HT03-RPL1C-Y

AERotex® - PES



Material	Colour	Size	Dimensions (W x H x D)	Weight	Equipment	Item No.
AERotex® - PES (Polyester)	■ Signal Red (RAL 3001)	L	approx. 36 x 47 x 26 cm	approx. 2 kg	Modular Bag Set Polyester, 5-colour	④ HT01-RPL1C-R
	■ Ultramarine Blue (RAL 5002)					⑤ HT01-RPL1C-B
	■ Black					⑥ HT01-RPL1C-BLK

Modular Bag Sets



Material	Colour	Item No.
AERotex® - PLAN (Tarpaulin)	■ Signal Red (RAL 3001)	① HT03-RPL1-C-SET-R
	■ Luminous Yellow (RAL 1026)	② HT03-RPL1-C-SET-B
	■ Ultramarine Blue (RAL 5002)	③ HT03-RPL1-C-SET-Y
AERotex® - PES (Polyester)	■ 5-colour	④ HT01-RPL1-C-SET