# AEROcase®

## Emergency Bags | Emergency Backpacks

#### AEROcase® - Pro

## We're offering more quality and manufacturing know-how!

AEROcase® - Pro<sup>1R</sup> und Pro<sup>EMS</sup> emergency bags and backpacks provide a safe transport of equipment and are reliable partners when it comes to first aid.

A compact and lightweight design allow an easy and comfortable transport which makes them the ideal companion for locations with long distances or restricted rooms.

Easy accessible outside pockets and well structured interiors allow quick access to the equipment on site while highly retroreflective reflective strips provide safety during day and night.


The emergency bags and backpacks are extremely durable and well equipped for high loads in everyday's use due to their high-quality materials and manufacturing process.

#### **Versatile field of application**

- Medical service
- SEG
- Ambulance
- Companies
- Special forces
- Schools
- Civil protection
- Private

#### AEROcase® - Pro Material details

Materials for highest requirements and durability



#### AEROtex® - PES

such as tear- and abrasion resistance.

AEROtex® - PES stands for high polyester 600 D (600 x 600). The nylon fibers are PVC-coated on the inner surface and are extremely robust in order to give the material all necessary properties for the daily use



#### AEROtex® - Plan

AEROtex® - Plan is a high quality and robust tarpaulin material.

It was designed for professional use in emergency medical service especially. Due to the both-sided PVC coating, the tarpaulin material is very resistant to loads and is furthermore waterproof. It is characterized by its scratch-, tearand abrasion restistance as well as being very hard-wearing. The surface can be cleaned with appropriate cleaners for tarpaulin material and appropriate disinfectants.

#### The quality in detail:

#### The improvements at a glance

we are happy to herby inform you about major improvements and design changes of our popular bags and backpacks from our **AEROcase® - Pro<sup>EMS</sup> Emergency Medical Service** series.

- Improved tarpaulin material: desinfectable, tear and abrasion resistant
- Improved polyester 600 D (600 x 600) and nylon (smooth, soft, durable)
- · High-quality stitching and product finish
- Improved stiff plastic handels in black and yellow and as well the soft-padded hook-and-loop handle
- Handles enforced and fastened with robust rivets
- 3 New comfortable extensions in loop shape!


- Soft-padded net in back area
- Slimmer fitted luminescent strips
- Quality with recognition factor







Only HT01/03-RPL1C-...







## AEROcase®

## **R** First Responder

### **Emergency Backpack**

#### AEROcase® - Pro 1 First Responder PL1C Emergency Backpack L

Versatile applications, high functionality and high-quality processing qualify the Pro<sup>1R</sup> PL1C for the perfect standard emergency backpack for professional use.

Flexibility and individual opportunities for organisation.

#### **Processing details**

- Eye-catching zippers for a better cognition
- Velcro straps (4 cm) for patches
- · Reflective strips all around
- Transparent pocket for a good visible placement of documents and information

#### **Professionally organized**

#### 2 front pockets

- Upper one for writing material
- The one below with plenty of space for ampoule cases I intubation cases etc.
- Additional zipper compartment for documents

#### 2 side pockets

• Space for sphygmomanometer, stethoscope, sharps container etc.

#### **Ergonomically formed handles**

- · Slip-proof comfort due to the reinforced core
- · Signal colors for a better cognition



Improved yellow standard handle
- for red/blue model -



Improved black standard handle
- for black model -



#### The models

Dimensions: 36 x 47 x 26 cm

Weight: 

2 kg | 3 kg

Volume: 

approx. 23 L













#### More equipment details

#### **① Clearly arranged interior concept**

The modular and structured interior concept ensures a secure placement of the equipment.

- Fold-out front
- Main pocket contains removable frame with extra seperating strip
- 5 module cases for the interior:
  - 3 in their corresponding divided main compartment
  - 2 can be placed inside the front pocket
- Elastic attachment straps (quality A)
- Lockable main pocket through cover flap with hook-and-loop fastener

#### **2** Comfortable carrying system

- Padded shoulder straps
- Additional quick release belts for fixation on stomach and upper body

#### **③ Protected all around**

- Dimensionally stable exterior walls
- Reinforced bottom protection of tarpaulin material
- · Robust plastic pyramid feets



Net in back area

### **Individual organisation**

- 5 modular inner bags in different sizes
- Tarpaulin 1 color (red, blue or yellow)
- Polyester assorted colors
- · Clear window and handles







Model: Polyester

