



## AERObase® – SEG XL of the EMS Series – New Product



### HT03-MPSEG1C-Y

(Luminous Yellow - RAL 1026)

Dimensions: approx. 37 x 54 x 17 cm

Weight: approx. 3.26 kg (incl. modular bags)

Volume: approx. 34 l

### Description:

Extremely durable SEG emergency backpack of the EMS series in "Extra Large" in a slim design. Made of AERObase® – PLAN (Tarpaulin), available in different colours as well as in black with a special matt finish. Developed for the continuous use in the emergency medical service and for MANV scenarios. Safe transport of emergency equipment with the highest possible carrying comfort and with large space. Delivery including 4 modular bags.

### Special features at a glance

- Interior with 4 modular bags with hook and loop fastening
- Inside of the lid with various compartments and loops, suitable for documentation such as medical emergency triage tags (Met-Tags), for small equipment and pens
- With flat front compartment, suitable for e.g. clipboard
- Robust, non-slip, disinfectable carrying handles
- Transparent window for inserting labels
- Extra strong and double riveted metal eyelets for back carrying straps (colour markings R/L alignment)
- Easy-grip, disinfectable zipper pullers in Luminous Yellow and optionally available in Grey HT10-ZP02C-GR/10 – SU 10 pcs.

### Further models



### HT03-MPSEGXL1C-R

Signal Red - RAL 3001



### HT03-MPSEGXL1C-B

Ultramarine Blue - RAL 5002



### HT03-MPSEGXL1C-BLK

Jet Black - RAL 9005, with a special matt finish

**Do you have questions or feedback?**  
**We are always there for you!**

Best regards,  
HUM Gesellschaft für Homecare und Medizintechnik mbH